

## Testing on the isokinetic dynamometer

Testing of joints - basic protocol	5.000,00
Testing of joint - trunk protocol	7.000,00
Testing in the closed kinetic chain - both extremity at the same time	5.000,00
Testing in the closed kinetic chain - both extremity separately	7.000,00
Adding a new modality to the basic protocol	2.000,00
Adding a new angle to the basic protocol	2.000,00
Adding a new regime to the basic protocol	2.000,00
Adding a new levels to the basic protocol	3.000,00
Adding a new joints to the basic protocol	3.000,00
Expert analysis and processing of data	3.000,00

## Training and rehabilitation on isokinetic dynamometer

Training and rehabilitation of joints/trunk - 8 trainings	45.000,00
Training and rehabilitation of joints/trunk - 12 trainings	60.000,00